


1. Introduction

A two dimensional HETeronuclear CORrelation experiment (HETCOR) is described here as ¹³C-¹H Correlation by Polarization Transfer. The experiment encodes the proton chemical shift information into ¹³C signals that are observed, and yields cross signals for all protons and ¹³C nuclei that are connected by a ¹³C-¹H coupling over one bond. The assignment of one member of a spin-coupled pair immediately leads to the assignment of the other. The HETCOR experiment described here uses a phase cycling for magnitude processing. This example demonstrates the basic principle of 2D NMR data acquisition and processing on Tecmag spectrometers.

2. Pulse sequence



Event Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Name:	Phrst	unblk	H90	t1/2	C180	tau/2	d2	H90'	C90	H90'	d3	ringdown	rx on	acq.	relax
Delay	1u	2u	H90	tau/2	C180	tau/2	d2	H90'	C90	H90'	d3	rd	ad	Acq. Time	Last Delay
F1_Ampl		H90 amp	H90 amp					H90 amp	H90 amp	H90 amp		F1 amp	F1 ...	F1 amp	
F1_PhMod												AS			
F1_Ph		X	X					phH2	phH2	phH2					
F1_Atten		H90 attn	H90 attn					H90 attn	H90 attn	H90 attn		F1 attn	F1 attn	F1 attn	
F1_TxGate															
F1_PhRst															
F1_UnBlank															
Acq															
Acq_phase														phrx	
RX_Blank															
RX_PhRst															
F2_Ampl						F2 amp			F2 amp						
F2_PhMod															
F2_Ph								phC1		phC2					
F2_Atten						F2 attn			F2 attn						
F2_TxGate															
F2_PhRst															
Delay_2D				de1:2				de1:2							

Acquisition	Frequency	Multi Rec.	Processing	Grad. Preemph.	Misc.	Sequence
H90	42.5u	H90'	=[(H90)[C90]]/2	ad	5u	F1 amp 90
tau/2	1u	C90	10.5u	Acq. Time	61.44m	H90 attn 14
C180	21u	d3	2.29m	Last Delay	3s	F1 attn 18
d2	3.45m	rd	25u	H90 amp	42.5	F2 amp 100

1D phase tables:

- (P1): 0.
- ph1(P2): 0, 2, 1, 3.
- ph0(P3): 0, 0, 0, 0,
2, 2, 2, 2.
- ph1(P4): 0, 0, 0, 0,
1, 1, 1, 1,
2, 2, 2, 2,
3, 3, 3, 3.
- phrx: (0, 2, 1, 3)₂,
(1, 3, 2, 0)₂,
(2, 0, 3, 1)₂,
(3, 1, 0, 2)₂.

(All tables are in 4 step mode.)

2D delay table:

- de1:2 (t1/2):
Auto, Every Pass.
Increment value: 200 μs,


Fig. 1. (a) The pulse sequence for the HETCOR experiment. (b) The sequence in the NTNMR sequence editor with phase tables set for magnitude processing.

3. Experiment

- Sample: 20 % ethyl crotonate in CDCl₃
- Spectrometer: 7 Tesla magnet with Tecmag HF3 Discovery
- Probe: Nalorac D300-5 OWB 5mm ¹H/¹³C Switchable probe
- ¹³C 90° pulse: 10.5 μs
- ¹H 90° pulse: 42.5 μs

3. Experiment (continued)

SW +/- (1D, ^{13}C) :	8.3 kHz (110.5 ppm)
SW 2D (+/-, ^1H):	1.25 kHz (4.15 ppm)
Dwell time (1D):	60 μs
Dwell_2D:	400 μs
Acq. points (and Points 1D):	1024
Points 2D:	256
Scans 1D:	32



Ethyl crotonate

Notes:

1. The sample is locked but not spun. The magnet is shimmed to ~ 0.8 Hz line width.
2. First, calibrate the 90° pulse widths of ^1H and ^{13}C using the nutation experiment (see note, "One Pulse Experiment and Pulse Calibration").
3. Set up the WALTZ sequence according to the note, " ^{13}C NMR Spectra with ^1H WALTZ Decoupling".
4. The center of pulses P2 and P4 should be aligned. Since $P2 > P4$, P2 have to split into 3 pulses. The delay of P2's middle pulse equals to P4, and the delay of both sides is $(P2 - P4)/2$. The middle pulse of P2 falls on the same event as P4 (and P5).

4. Results

Data processing:

In the NDFT window, click the "1D" tab and select: "Use 1D settings", "Fourier Transform", "Sine Bell, SB Shift: 90, SB Skew: 5, SB width: 1024"; click the "2D" tab and select: "Use 2D settings", "Zero fill, 1", "Sine Bell, SB Shift: 90, SB Skew: 5, SB width: 512"; "Fourier Transform", "Transpose: complex", click "Do it". □

The data processing procedure is automated in an NMRScript. To use the script, click "Scripts|Processing scripts|Process 2D (Magnitude)".




Fig. 2. Magnitude plot of the HETCOR spectrum of ethyl crotonate.

5. References

1. R. Freeman, G. A. Morris, *J. Chem. Soc. Chem. Commun.*, **1978**, 684-686.
2. S. Braun, H.-O. Kalinowski, S. Berger, "150 and More Basic NMR Experiments", Wiley-VCH, 1999, 375-377.