

1. Introduction

A **Two Pulse Phase Modulation (TPPM)**<sup>[1]</sup> proton decoupling scheme is often used to improve proton decoupling in double resonance magic angle spinning (MAS) experiments. Figure 1 shows a TPPM scheme following a cross-polarization from proton to a rare spin. TPPM is conveniently implemented with an asynchronous event incorporating a two-event sequence.

2. Pulse sequence




Fig. 1a: CP-TPPM pulse sequence. b: Actual pulse sequence in the NMR sequence editor.

### 3. Experiment

Sample:  $u\text{-}^{13}\text{C}$ ,  $^{15}\text{N}$ -histidine

Rotor speed: 7.5 kHz

Proton  $90^\circ$  pulse width: 4.5  $\mu\text{s}$

Mixing  $rf$  field: 55.5 kHz

Proton decoupling field: 62.5 kHz

TPPM flip angle =  $150^\circ$  ( $\tau = 8.3 \mu\text{s}$ );  $\phi/2 = 12.5^\circ$

Number of scans: 4




Fig. 2.  $^{13}\text{C}$  CPMAS spectra of  $u\text{-}^{13}\text{C}$ ,  $^{15}\text{N}$ -histidine with (a)  $cw$  decoupling and (b) TPPM decoupling. Note improvement in linewidth (75 Hz to 145 Hz) and corresponding increase in signal to noise ratio.

### 4. Reference

1. A.E, Bennett, C.M. Rienstra, M. Auger, K.V. Lakshmi, and R.G. Griffin, *J. Chem. Phys.* 103, 6951, 1995.